

**Role of mycobacterial protein, PPD in modulating TLR
induced T-cell response**

A Dissertation Thesis

Submitted to

Nirma University

in Partial Fulfilment of requirement for the

Degree of

Master of Science

in

Biotechnology/Biochemistry



Submitted by

Jagdish Madariya (10MBT004)

Prerak Vyas (10MBC020)

Under the guidance of

Prof. Sarat K. Dalai

ACKNOWLEDGEMENT

We express our deep sense of gratitude towards the **Almighty God** for providing us the strength, energy and inspiration to work on our project. We would also like to thank our family member for constant support for our project.

We convey our special thanks to **Prof. Sarat Dalai** whose vision, motivation and constant support helped us carry out our project and for all the painstaking efforts and deep insights into the problems which helped us troubleshoot them at any hour of the day; improving the quality of our work at all stages. No words can express our sincere and deep sense of reverence for him. We are extremely indebted to him for the scientific attitude he has installed in us which will definitely stand in all future endeavours and it was because of him that we were able to learn so much in this short period. We thank him from the bottom of our soul to teach us the lessons to be good human beings and leading our lives with sense of patience, satisfaction and enthusiasm.

We would also like to thank other faculty members Dr. Shalini Rajkumar, Dr. Ameer Nair, Dr. Sonal Bakshi, Dr. Nasreen Munshi, Dr. Shriram Seshadri, Dr. Mili Das and Dr. Vijay kothari for their extensive support in our studies and guiding us for such opportunity.

We are also very grateful to Ph.D. scholars especially Mr. Hardik Patel for giving us constant guidance regarding our project and lab work. We extend our thanks to Mr. Prasant, Mr. Rahul, Mrs. Kaveri, Ms. Harshita, Mrs. Disha, Mr. Maharshi, Mr. Mahendrapal, Mr. Palakl, Mr. Arpan, Ms. Purvi, Ms. Arpita and Ms. Suhani for their support. We would also like to thank Rimanshu for his help for our project. We express our sincere thanks to Mr. Sachin Prajapati, Mr. Bharat Anand and Mr. Rajendra for providing with all necessary equipments, glasswares and chemicals at any time the day during our work.

We are also thankful to **Dr. Hemangini Vora**, Assistant Scientist, Department of Immunohistochemistry, Gujarat Cancer Research Institute (GCRI), Ahmedabad, for providing access to Flow Cytometry for T- cell assay. We would also like to thank Birva, Research scholar at GCRI for constant support in taking reading on Flow Cytometer.

We extend our special thanks to all our friends in our dissertation group-, Arva Sharaf, Krunal Chodvadiya, Parini Surti, Depika Mandaliya, Manoj Patidar and Naveen Yadav who worked in coordination with us and provided appreciable moral support throughout dissertation period and helped us in solving various problems that we

encountered during our project. We would also like to extend our sincere thanks to our **beloved parents** for blessing us on the occasions when we needed them the most.

In the end our deepest gratitude to all those who have not been mentioned by name but nevertheless have been of invaluable help in their inscrutable ways.

Jagdish Madariya and Prerak Vyas

Date:

Place: Ahmedabad

(Jagdish Madariya) (Prerak Vyas)

CONTENTS

1. ABBREVIATIONS

2. LIST OF TABLES

3. LIST OF FIGURE

4. INTRODUCTION

5. MATERIALS AND METHODS

6. RESULTS AND DISCUSSION

7. SUMMARY

8. REFERENCES

9. APPENDIX

ABBREVIATIONS

Ag: Antigen

AP-1: Activator Protein 1

APC: Antigen Presenting Cell

Bcl-2: B cell Leukemia 2

BM: Bone marrow

CD: Cluster of Differentiation

CM: Central memory

DC: Dendritic Cell

DEPC: Diethyl pyrocarbonate

EDTA: Ethylene Diamine Tetra acetate

ELISA: Enzyme Linked Immuno Sorbant Assay

EM: Effector memory

EtBr: Ethidium Bromide

FACS: Fluorescence associated cell sorter

IFN- γ : Interferon γ

IL-15: Interleukin 15

IL-15R: Interleukin 15 Receptor

IL-7: Interleukin 7

IL-7R: Interleukin 7 Receptor

JAK: Janus kinase

LLCM: Long Lived Central Memory

LPS: Lipopolysaccharide

LTA: Lipoteichoic acids

MAPK: Mitogen Activated Protein Kinase

MHC: Major Histocompatibility Complex
mIL- 15: Mouse(Mus musculus) Interleukin 15

mIL -7: Mouse(Mus musculus) Interleukin 7

MOPS: 3-[N-morphino]- propanesulfonicacid

MyD88: Myeloid Differentiation Factor 88

NFKB: Nuclear Factor Kappa B

NK cells: Natural Killer cells

NKT cells: Natural Killer like T cells

PAMP: Pathogen Associated Molecular Pattern

PBS: Phosphate Buffer Saline

PCR: Polymerase Chain Reaction

PGN: Peptidoglycan

PI3K: Phosphatidylinositol-3 kinase

PKB: Phosphokinase B

Poly I: C: Poly inosinic: cytidinic acid

PRR: Pattern Recognition Receptor

RT: Reverse Transcription

RT-PCR: Reverse Transcription Polymerase Chain Reaction

STAT: Signal Transducer and Activator of Transcription

TAE buffer: Tris Acetate EDTA buffer

TCR: T cell Receptor

TH cells: T helper cells

TLR: Toll like receptors

TIR: Toll-like IL-I Receptor

Tm: Melting Temperature

TNF: Tumor Necrosis Factor

TRAF: TNF Receptor Associated Factor

TRAF: TNF Receptor-associated factor

TU: Tubercullin Unit

LIST OF TABLES

Table no.	Title	Page number
1	Toll like receptors and their ligand(TLR ligand)	06
2	Immunization protocol for first objective	14
3	Immunization protocol for second objective	15
4	β -actin primer detail	21
5	IL-15 primer detail	21
6	Composition of Reaction Mixture for RT	24
7	Composition of Reaction Mixture for PCR	24
8	PCR profile for β -actin	26
9	PCR profile for IL-15	26
10	% of dividing CD4 ⁺ T-Cells with respect to total dividing cell (TDC) and Total cell (dividing+ non dividing cells)	39
11	% of dividing CD8 ⁺ T-cells with respect to total dividing cell (TDC) and Total cell (dividing+ non dividing cells)	41
12	IL-15 ELISA	45
13	% of dividing CD8 ⁺ T-Cells with respect to total dividing cell (TDC) and Total cell (dividing+ non dividing cells)	53
14	% of dividing CD4 ⁺ T-Cells with respect to total dividing cell (TDC) and Total cell (dividing+ non dividing cells)	55
15	IL-15 ELISA	59
16	Absorbance of RNA at 260nm, 280nm and A_{260}/A_{280}	61

LIST OF FIGURES

Figure No.	Title of the figures	Page No.
1	Kinetics of Innate and adaptive immunity	2
2	Stimulation of adaptive immunity by innate immune response	3
3	Development of the cells of immune system	3
4	Phases of T-cell responses	4
5	Cytokines control T-cell differentiation at different phases	5
6	TLRs and TLR ligands	6
7	Pathogen recognition of TLRs and processing and presentation of antigen to T-cell	6
8	cytokine and chemokines induction by CFA	9
9	Mechanism of action of CFSE	16
10	CFSE dilution with cell division	16
11	Scheme of analysis of Flow cytometry data	30
12	T-cell responses in mice immunized with OVA alone	31
13	T-cell responses in mice immunized with OVA in presence of PPD	32
14	T-cell responses in mice immunized with OVA and LPS	33
15	T-cell responses in mice immunized with OVA,LPS and PPD	34
16	T-cell responses in mice immunized with OVA and zymosan	35
17	T-cell responses in mice immunized with OVA,zymosan and PPD	36
18	T-cell responses in mice immunized with OVA, LPS and zymosan	37
19	T-cell responses in mice immunized with OVA, LPS and zymosan and PPD	38
20	Bar diagram of % of dividing Ova specific CD4 ⁺ T-cell	40
21	Bar diagram of % of dividing Ova specific CD4 ⁺ T-cell	41
22	Bar diagram of % of dividing Ova specific CD8 ⁺ T-cell	43
23	Bar diagram of % of dividing Ova specific CD8 ⁺ T-cell	43
24	Line graph of % of dividing Ova specific CD4 ⁺ T-cell	44

25	Line graph of % of dividing Ova specific CD4 ⁺ T-cell	44
26	IL-15 ELISA of second objective	46
27	T-cell responses in mice immunized with OVA alone.	47
28	T-cell responses in mice immunized with HEL alone	48
29	T-cell responses in mice immunized with OVA and HEL	49
20	T-cell responses in mice immunized with OVA, LPS and zymosan	50
31	T-cell responses in mice immunized with OVA, LPS and zymosan and HEL	51
32	T-cell responses in mice immunized with HEL, LPS and zymosan	52
33	Bar diagram of % of dividing Ova specific CD8 ⁺ T-cell	54
34	Bar diagram of % of dividing Ova specific CD4 ⁺ T-cell	56
35	Number of division CD4 ⁺ T-cells has undergone	57
36	Number of division CD4 ⁺ T-cells has undergone	57
37	% of dividing CD4 ⁺ T-cells with increase in Ova challenge	58
38	IL-15 ELISA of second objective	59
39	Electrophoresis of isolated RNA sample	60
40	electrophoresis of cDNA transcribed from isolated RNA	61
41	electrophoresis PCR product for β -actin	62
42	RNA isolation from splenic dendritic cells stimulated with TLR ligands	63
43	Polymerase Chain Reaction of DC cells stimulated with TLR ligands (<i>In vitro</i>)	63

INTRODUCTION

1.1 Innate and adaptive immunity

1.2 Development of the cells of immune system

1.3 T-cell activation and differentiation:

1.4 1.4 TLR ligand and their receptor

1.5 Adjuvants (IFA and CFA)

1.6 Purified protein derivatives (PPD)

1.7 TLR ligands used in our study

a) LPS

b) Zymosan

INTRODUCTION

1.1

Host defence against invading microbial pathogens is elicited by the immune system, which consists of two components: innate immunity and adaptive immunity. The innate immunity provides the first line of defence against infection. It consists of cellular and biochemical defence mechanisms that are in place even before infection and are poised to respond rapidly to infections.

In adaptive immunity, B and T lymphocytes utilize antigen receptors such as immunoglobulins and T cell receptors to recognize non-self. Adaptive immune response against an antigen is generated within five or six days after the initial exposure




Figure1: Kinetics of Innate and adaptive immunity

to that antigen. It responds to the infection with a high degree of specificity as well as the remarkable property of “memory”, that provide long-lasting immunity against the re-infection. Because the adaptive immune system takes longer to be activated pathogens get favourable condition for their rapid growth. During that period the innate immune system takes responsibility by acting immediately to encounter the infections. Thus innate immunity mainly performs two function. One is the initial response to microbes (antigen) and second is stimulation of adaptive immunity by processing and presenting the antigen to T-cell. This processing and presentation is done by APC (antigen presenting cells like macrophages and dendritic cells) The innate and adaptive

immune system is linked through direct cell contact and through interactions involving chemical mediators such as cytokines and chemokines.




Figure 2:- Stimulation of adaptive immunity by innate immune response

(Cellular and molecular immunology by Abul K. Abbas 6th edition)

1.2 Development of the cells of Immune system:

All cells of the immune system have their origin in the bone marrow and they include myeloid (neutrophils, basophils, eosinophils, macrophages and dendritic cells) and lymphoid (B lymphocyte, T lymphocyte and Natural Killer) cells, which differentiate along distinct pathways




Figure 3: Development of the cells of immune system

Bone marrow precursor T-cells enter the thymus through blood circulation and undergo the process of maturation and become a naive T-cell. After the activation these naive T-cells differentiate into antigen specific effector T-cells

1.3 T-cell activation and differentiation:

Bone marrow precursor T-cells migrate to thymus and undergo the process of maturation to become naive T-cells. After maturation these naive cells moves to secondary lymphoid organ, where they encounter antigen presented by mature dendritic cells on class I and II MHC molecules and thus become activated. This results in the expansion of the antigen specific lymphocytes the differentiation of these cells into effector and memory T cells




Figure 4: Phases of T-cell responses

When APC present antigen to T-cell by MHC molecule, they also provide a complex set of co-stimulatory signal to T-cell and helps them to undergo vigorous clonal expansion and acquire effector function. In adaptive immunity, B and T lymphocytes utilize antigen receptors such as immunoglobulins and T cell receptors to recognize non-self. Adaptive immunity takes five to six days to generate response against antigen. But after activation they form memory cells which provide long lasting immunity against the re-infection.

Effector CD4⁺ T-cell thus produced respond to antigen by producing cytokines that have several actions like activation of macrophages and B lymphocytes, and CD8⁺

CTLs respond by killing other cells. These cytokines control events of T-cell activation at various stages of T-cell responses.

There are three cytokines- IL-2, IL-7 and IL-15 which have been shown to have integral roles in CD8+ T-cell memory generation and maintenance. Although their functions might be overlapping in some cases, their activities on T-cells can occur either simultaneously or at temporally distinct points during an immune response




Figure 5: Cytokines control T-cell differentiation at different phases.

The antigen in our system is recognised by the cells of innate immunity like APC and dendritic cells. These cell recognise the specific pattern on pathogens i.e. PAMPs (pathogen associated molecular patterns) with the help of receptor molecules present on their membrane called Toll like receptors(TLRs)

1.4 Toll Like Receptors (TLRs)

TLRs are membrane-bound molecules that recognize microbial components on the surface or within extracellular compartments of antigen presenting cells (APCs). TLRs are pattern recognition molecules, which recognise specific pattern on pathogen and activate innate immune system, which in turn stimulate adaptive immune system.

These TLRs binds to TLR ligand like lipopolysaccharide (LPS) from Gram-negative bacteria. LPS is a TLR ligand which binds to TLR4. Till now at least 11 mammalian TLRs have been identified among which TLR1-9 are conserved between the human and mouse.




Figure 6:TLRs and TLR ligands

TLR Ligands are the microbial components which bind with different TLRs and activate specific TLR pathways. Stimulation of TLRs by microbial components triggers expression of several genes that are involved in immune responses, like genes for inflammatory cytokines and co-stimulatory molecules

Table 1: Toll like receptors (TLRs) and their ligand(TLR ligand)

S.no	Toll like receptors (TLRs)	TLR ligand recognised by TLRs
1	TLR2+TLR6	diacyl lipopeptide
2	TLR2+TLR1	triacyl lipopeptide
3	TLR4	Gram –ve LPS,fungal mannans, parasitic phospholipids
4	TLR5	Bacterial flagellin
5	TLR9	CpG DNA
6	TLR3	Virsl ds RNA




Figure 7:Pathogen recognition of TLRs and processing and presentation of antigen to T-cell

The molecular mechanism lying behind the expression of various gene is the dimerization of TLRs. TLR2 is shown to form a heterophilic dimer with TLR1 or TLR6, but in other cases TLRs are believed to form homodimers. Dimerization of TLRs triggers activation of signaling pathways, which originate from a cytoplasmic TIR domain.

TLRs leads to expression of many gene that give signals to other cells of immune system to respond the antigen in a better way. This is the basis on which adjuvants works and helps in enhancing the immune response against pathogens.

1.5 Adjuvants:-

The word adjuvant comes from the Latin word adjuvare, which means to help or to enhance. The role of innate immunity in stimulating adaptive immune responses is the basis of the action of adjuvants, which are compound that when administered with antigen potentiates adaptive immune response, both humoral & T-cell mediated response, against co-administered antigen. The main function of adjuvants is to induce inflammatory response that will help in triggering adaptive immune response. In absence of adjuvant even higher dose of antigen cannot induce productive immune response, which can be achieved with low amount of antigen if adjuvant is used.

Adjuvants are useful in experimental immunology and in clinical vaccines. Many adjuvants in experimental use are microbial products, such as **killed mycobacteria** and **LPS**, that engage TLRs and elicit strong innate immune responses at the site of antigen entry.

Alum is one of the few adjuvants approved for human use in most of the countries worldwide. Alum induce a good antibody (Th2) response but it has a little capacity to stimulate cellular(Th1) immune responAse which is so important for protection against many pathogens. Alum is thought have a role in neurodegenerative diseases such as Alzheimer's disease. Pure recombinant or synthetic antigens used in modern day vaccines are generally far less immunogenic than older style live or killed whole organism vaccines. So there is a need for safer and more effective adjuvants suitable for human use.

Since last few decades, adjuvants like IFA (Incomplete Freund's adjuvant) and CFA (Complete Freund's adjuvant) are in use for experimental work and preclinical vaccine preparation to enhance immune response. Its mode of action, however, is still not very clear.

IFA is essentially a paraffin oil containing mannide mono-oleate as a surfactant.. When IFA is mixed with aqueous solution of antigen it forms water in oil emulsion where antigen is in water phase. This water in oil emulsion slowly releases antigen and help to induce immune response effectively (Alfonso Billiau et.al 2001). CFA contain an additional factor to oil i.e. heat killed *Mycobacterium tuberculosis*. It is more potent than IFA. It is used in primary immunization to potentiate adaptive immune response by triggering innate immunity

IFA stimulates innate immunity, as evident by transiently increased resistance to bacterial infection (Castro et. al., 1993).It also induces the expression of cytokines, predominantly tumor necrosis factor α (TNF- α) in regional lymph nodes (Mussener et. al., 1995). Besides inducing the above responses, CFA containing mycobacteria attracts macrophages and other cells to site of injection and enhances immune response.

Various components of heat killed mycobacteria interact with various TLR which activates innate immune system and attracts macrophages and other cells to the site of infection and enhances immune response

Active components of CFA:

- Heat shock protein (HSPs), the intracellular protein found in both prokaryote and eukaryote, bind to TLR4, resulting in an LPS like effect on mononuclear phagocytes (Kol et. al., 2000; Ohashi et. al., 2000).
- Glycolipids typical for mycobacteria trehalose dimycolate (TDM) and lipoarabinomannan (LAM). LAM is the agonist of the TLR2 (Golenbock et. al., 2001).
- CpG oligodeoxynucleotides, present in heat killed mycobacteria have been shown to strongly stimulate innate immune mechanisms, including production of cytokines by mononuclear phagocytosis (MPCs), maturation, and activation of APCs (Weiner, 2000).
- Muramyl-di peptide (MDP) is also immuno-modulator component of CFA which is the building block of peptidoglycan component of the bacterial cell wall. MDP was found

in expansion of granulocyte-macrophage progenitors in spleen, and proliferation of multi potential stem cells in bone marrow (Galelli et. al., 1983).

Mechanism of Immuno-potentiality by CFA:-

- **Enhancement of antigen uptake by APCs:** Water in oil emulsion of antigen causes slow delivery of antigen to immune system, hence increased persistence of antigen in system.
- **Complex set of signal:** The mycobacteria in CFA are currently assumed to be recognized by PAMP (pathogen-associated molecular pattern) receptors on various immunocompetent cells and thus to provide the stimulus for these cells to release mediators and express membrane receptors (collectively designated as “danger” signals) that will lead to Th1-type skewing of ongoing immune responses
- **Cytokine and chemokine induction:** Cytokines observed to be induced in the early phases following exposure to CFA (or mycobacteria) are TNF- α , IL-12, IL-6, IFN- γ , and several chemokines. Mycobacterial components are known to target MPCs and DCs (involving TLRs) and to induce production of monokines, in particular IL-12 and TNF- α (Mussener et. al., 1995).

CFA is a more potent adjuvant because it contains heat killed *mycobacterium*. The antigens present on *mycobacterium* are mainly responsible for the generation of enhanced immune response. PPD (purified protein derivatives) is one of the major mycobacterial protein antigen. These antigens boost the response by providing complex set of signals to immune cells.




Figure 8: cytokine and chemokines induction by CFA

1.6 PPD (Purified protein derivatives):

Tuberculin Purified protein derivatives (PPD) is an antigenic preparation derived from *Mycobacterium tuberculosis*, the causative agent of tuberculosis. It has had a widespread clinical and epidemiologic, and investigative use. It is used as an in vivo diagnostic aid for detection of mycobacterium tuberculosis infection. It is basically used in the measurement of T cell reactivity in persons exposed to TB or vaccinated with BCG.

PPD is prepared by precipitation of proteins from heated *Mycobacterium tuberculosis* cultures and composed of soluble protein antigens.

PPD contains multiple antigenic component. One such component, a small protein or a large peptide designated as antigen 7, was identified as major constituent of PPD. D-Arabinose and D-mannan is major polysaccharide constituent of PPD. (Yu Ma and Thomas M. Daniel JSTOR 1983).

PPD stimulate innate immunity (Jenefer M. Blackwell et al. 2004); elicit strong Th1-cell response which can be identified by production of IFN- γ , IL-2, and TNF- α , TNF- β . PPD also act as immunogen carrier (De Silva BS 1999) and helps in inducing effective adaptive immune response. It has potential in vaccine preparation as an adjuvant. (De Silva et al., 1999). The effect of PPD in modulating TLR induced T-cell response is not yet clear.

In the proposed study we want to test whether enhanced antigen specific immune response in presence of PPD is due to modulation of TLR induced response of adjuvants or due to modulation of innate immunity by TLR ligands.

1.7 TLR Ligands used in our study:

We have used following two TLR ligand in different combination with antigen.

- **Lipopolysaccharide (LPS)**

LPS is the major constituent of cell wall of Gram negative bacteria. It is a potent immunostimulator. It is a TLR4 ligand recognised by TLR4 receptor.

LPS stimulates APCs to release various cytokines e.g., IFN, TNF- α , IL-12 and also elicit both CD4⁺ and CD8⁺ T-cell response (David et al., 1997). HSP60 elicits potent inflammatory response in innate immune cells. HSP60 signalling mediated by TLR4 (LPS) release inflammatory mediators such as IL-6, IL-15, TNF- α (Ohashi et al., 2000). LPS induces upregulated expression of both IL-15 and IL-15R α by splenic

dendritic cells and stimulation of antigen specific T cell proliferation (Mattei et. al., 2001).

- **Zymosan**

Zymosan is a yeast cell wall particle. It is a TLR2 ligand recognised by TLR2 receptor. It has potential to activate macrophages, monocytes, and leukocytes resulting in the stimulated secretion of inflammatory products including TNF- α , IL-8, hydrogen peroxide, and arachidonic acid. Zymosan also induce NF- κ B activation(Sato et. al., 2003). It is also able to induce maturation of dendritic cells *in vitro* and to stimulate the production of IL-12, provide a link between innate and adaptive immune response.

We want to mimic CFA by using by using LPS (which represent TLR4), zymosan (which represent TLR2 of mycobacteria) and PPD which is a major mycobacterial protein antigen.

When two or more antigens are injected together they can respond in either way, they enhance the immune response against each other or they dampen the response against each other. But it is well known that in the case of multivalent vaccine like BCG two or more antigen mixed together enhances the immune response against each other.

So Further we want to test the idea whether one antigen from a pathogen can influence the immunogenicity of another antigen.To test this idea we have taken two protein antigen chicken Ovalbumin and hen egg white lysozyme (HEL) which are the from same source

MATERIALS AND METHODS

2.1. MICE

2.2. TLR LIGANDS

- **LPS**
- **ZYMOSAN**

2.3. ANTIGENS

- **OVALBUMIN**
- **PPD**
- **HEN EGG WHITE LYSOZYME (HEL)**

2.4. IMMUNIZATION PROTOCOL

2.4.1 PROTOCOL OF FIRST OBJECTIVE

2.4.2 PROTOCOL OF SECOND OBJECTIVE

2.5. T-CELL PROLIFERATION ASSAY

2.6. RT-PCR

2.6.1 BUFFERS AND REAGENTS

2.6.2 FOR CFSE DILUTION ASSAY

2.6.3 FOR RNA ISOLATION

2.6.4 FOR CDNA SYNTHESIS

2.6.5 FOR PCR AMPLIFICATION

2.6.6 FOR AGAROSE GEL ELECTROPHORESIS

2.6.7 FOR *IN VITRO* CULTURE

2.7. PRIMERS

2.7.1 MUS MUSCULUS ACTIN, BETA (ACTB), MRNA

2.7.2 MUS MUSCULUS INTERLEUKIN 15 (IL15), TRANSCRIPT VARIANT 1, MRNA

2.8. PROTOCOLS

2.8.1 RNA ISOLATION

2.8.2 CDNA SYNTHESIS

2.8.3 PCR AMPLIFICATION

2.8.4 IL-15 INDUCTION

2.1 Mice

Male BALB/C mice (Age: 4-5 week old) were obtained from Zydus pharmaceuticals, Ahmedabad. These animals were housed in Animal house facility, Institute of Pharmacy, NIRMA University. Food and water were available ad libitum. All experiments performed in accordance with protocol approved by Animal care and use committee of NIRMA University.

Immunizations of Mice

Incomplete Freund's adjuvant (IFA) [Banglore Genei, Cat# FIA]

IFA was used for water in oil emulsion at 1:1 ratio.

2.2 TLR Ligands:

1. Lipopolysaccharide (LPS): (Invivogen tlr1-pelps:-LPS-EB)

Stock solution: 5 mg of LPS was dissolved in 1 ml of sterile 1X PBS.

(Stock concentration: 5 mg/ml)

Dose given: 10 µg/ mouse

2. Zymosan (Invivogen, cat. #tlr1-zyn):

5 mg of Zymosan was weighed and dissolved into 5 ml of PBS.

(Stock concentration: 1 mg/ml).

Dose given: 25 µg/ mouse

2.3 Antigen

• Ova Protein:

Chicken Ovalbumin (Merck, cat.# 32467 1gmcn)

Stock solution: 40 mg of Ova Protein was dissolved in 1.0 ml Sterile PBS.

(Stock concentration:40mg/ml).

Dose given: 200 µg/ mouse

- **PPD (Purified Protein Derivative):**

PPD was purchased from Span Diagnostics, Surat, Gujarat.

Stock Concentration: 10TU/0.1 ml

Dose given: 0.2 µg/ mouse

- **Hen Egg white Lysozyme (HEL):** Merck (4403-1GMCN)

40 mg of HEL was dissolved in 1ml sterile PBS.

Stock concentration:40mg/ml

Dose given: 100 µg/ mouse

2.4 Immunization Protocol

2.4.1 Objective 1:-Role of PPD in modulation of TLR induced T-cell responses.

Following combinations of TLR ligands with antigen and/or PPD were injected to

Mice intraperitoneally

Table 2: Immunization protocol for first objective

S.no.	Combination of TLR ligand	Serum collection for IL-15 ELISA		Number of mice sacrificed for T-cell assay	
		24 hr	72 hr	10 th day	10 th day repeat
1	OVA+IFA	3	3	3	3
2	PPD+OVA+IFA	3	3	3	3
3	OVA+LPS+IFA	3	3	3	3
4	PPD+OVA+LPS+IFA	3	3	3	3
5	OVA+ZYMOSAN+IFA	3	3	3	3
6	PPD+OVA+ZYMOSAN+IFA	3	3	3	3
7	OVA+LPS+ZYMOSAN+IFA	3	3	3	3
8	PPD+OVA+LPS+ZYMOSAN+IFA	3	3	3	3

2.4.2 Objective 2: whether one antigen from a pathogen can influence the immunogenicity of another antigen of same origin

To test this idea we have taken two antigen ovalbumin and Hen egg white lysozyme (HEL), which are having same origin and injected these antigens in different combination with TLR ligands

Table 3: Immunization protocol for second objective

S.no.	Combination of TLR ligand	Serum collection for IL-15 ELISA at 72 hours	Number of mice sacrificed for T-cell assay
1	OVA+ IFA	3	3
2	HEL + IFA	3	3
3	HEL+OVA+IFA	3	3
4	OVA+LPS+ ZYMOSAN +IFA	3	3
5	HEL +LPS+ZYMOSAN+IFA	3	3
6	HEL +OVA+LPS+ZYMOSAN+IFA	3	3

2.5 T-cell Proliferation Assay:

The measurement of T cell responses is an important parameter to determine status of immune system. Usually T cell assays are done by measuring:

1. Cell proliferation
2. Synthesis of Intracellular Cytokine or secreted cytokine

Cell dilution assay can be done by using either thymidine assay or CFSE dilution assay. The classical method of measuring T cell proliferation has been done via uptake of [3H] thymidine during the final hrs of 3-5 day culture. But this method demands handling of radioisotopes which is not very safe. Moreover the interpretation of results is very subjective. The difficulties and limitations of this method have been overcome by CFSE dilution assay

Carboxyfluorescein succinimidyl ester (CFSE) is an effective and popular means to monitor lymphocyte division. CFSE covalently labels long-lived intracellular molecules with the fluorescent dye, carboxyfluorescein. Thus, when a CFSE-labelled cell divides, its progeny gets half the number of carboxyfluorescein-tagged molecules and thus each cell division can be assessed by measuring the corresponding decrease in cell fluorescence via Flow cytometry.

CFDA-SE is carboxy-fluorescein-diacetate-succinimidyl-ester is a lipophilic non-fluorescent compound. It diffuses into the cells passively. The acetate groups are removed by esterase activity inside the cell yielding highly fluorescent compound, which binds random secondary groups of all cellular proteins and ultimately label the cells




Figure 9:-Mechanism of action of CFSE

The capacity of CFSE to label lymphocyte populations with a high fluorescent intensity of exceptionally low variance, coupled with its low cell toxicity, make it an ideal dye to measure cell division. Since it is a fluorescein-based dye it is also compatible with a broad range of other fluorochromes making it applicable to multi-color flow cytometry.

This method is based on analysis of CFSE dilution upon *in vitro* recall stimulation with antigen. Analysis of proliferation was combined with the use of monoclonal antibodies directed against the lymphocyte surface markers like CD4, CD8, CD3 etc. This assay helps analyzing proliferation of cells as they divide in response to antigenic or mitogenic challenge and their fate division.




Figure 10: CFSE dilution with cell division

Upon antigenic or mitogenic challenge, cells start to proliferate and by this they lose fluorescence in step-wise manner with each cell division. Reduction in fluorescence intensity can be quantified by flow-cytometry. Thus, there are many advantages to use CFSE cell proliferation including determination of proliferation of specific cell population.

Labelling cells with CFSE:

- Lymph node obtained from mice was minced in sterile 1X PBS and 1% glucose solution. The cells obtained were pellet down by spinning at 1000 rpm for 10 mins. Cells in pellet were treated with 1 ml RBC lysis solution and kept for 3 min. (Optional step if lymph node is used). Repeat the washing step using 1X PBS.
- Cells were pellet down and then dissolved in chilled complete RPMI and taken for counting under microscope. Based on the counting definite volume of cell was taken in a way to stain 3×10^6 cells with $5 \mu\text{M}$ CFSE.
- To stain the cells uniformly dye is not added directly to cells as dye being denser will not stain cells uniformly. For this reason 3×10^6 cells were made up to certain volume in one tube and 5mM CFSE dye made in another tube with prewarmed complete RPMI. Both the tubes were mixed together in a way that working concentration of dye is $5 \mu\text{M}$ (CFSE dye is photosensitive so avoid intense light during CFSE staining).
- Cells were incubated at 37°C and 5% CO_2 for 10 min during CFSE staining in CO_2 incubator.
- The staining is quenched by addition of 2 volumes of chilled complete RPMI 1640 to the stained cells. The cells were pellet down at 1000 rpm for 10 min. Pellet is resuspended in the complete RPMI and spin as before, repeating washing step twice.
- Finally cells in pellet were resuspended in media in a way that $100 \mu\text{l}$ of 3×10^5 cells are dispensed to the each well of culture plate. $100 \mu\text{l}$ of different concentrations of antigen (OVA or PHA) are added to 96 wells culture plate. The culture plate is incubated at 37°C and 5% CO_2 for 72 hours in CO_2 incubator.
- After 72 hours cells were washed from the wells with 1% FACS buffer and collected in eppendorfs and volume was made upto 1 ml.

Cells were pellet down by spinning at 1200 rpm for 7 mins and resuspended in minimum volume. Then $10 \mu\text{l}$ of 1:50 diluted normal mouse serum was added and

mixed by tapping the tubes and by pipetting. Cells were kept for 10 mins on ice. After which 10µl 1:75 times diluted antibodies were added, mixed well by pipetting and tapping the tubes and kept on ice for 30 mins.

Note: Ensure that the serum and antibodies are bound well, for this reason pipetting and tapping the tubes are crucial steps that would give the better separation in the results.

- Dilutions are made in such a way that serum is eventually diluted 1:100 times and antibodies are diluted 1:150 times in the suspension of cells.

Note: Keep eppendorf in dark while antibody staining

- Volume was made upto 1 ml with 1% FACS buffer and centrifuged at 1200 rpm for 7 mins. Pellet obtained was washed once again with 1% FACS spin same as above.
- Finally the pellet was resuspended in 500µl of 1% FACS. 50µl of 1% Formaldehyde was added for storage up to 24 hours at 4⁰C. Samples were acquired using flow Cytometer (BD Canto) at GCRI, Ahmedabad. The acquired results were analyzed using the software FlowJo (7.6.5).

2.6 Buffers and Reagents required:

2.6.1 For CFSE dilution assay:

1. Lymph node or spleen
2. Frosted slide for preparing single cell suspension
3. Sterile dissection tools and petriplate
4. Sterile complete RPMI media (with 10% FCS)
5. Sterile RPMI washing media (with 1% FCS)
6. CO₂ incubator
7. Swinging rotor clinical centrifuge
8. CFSE [Carboxyfluorescein succinimidyl ester (cat# no. 21888-25MG-F, Sigma)]
9. Phytohaemagglutinin (PHA)
10. Flow Cytometer (BD FACS Canto II at GCRI, Ahmedabad)
11. Flowjo software (version 7.6.5) to analyse result

12. 96 Well U-bottom culture plate.

13. Fluorescently labelled antibodies

- APC Rat Anti-Mouse CD4 (BD Bioscience, cat# 553932)
- PE Rat Anti-Mouse CD8a (BD Bioscience, cat# 554714)
- PerCP Hamster Anti-Mouse CD3 (BD Bioscience, cat# 553975)

2.6.2 *In vitro* culture

1. Phosphate buffer saline, PH 7.4 (PBS)

2. Complete RPMI 1640 media (10% FCS) (HI GlutaXLtm RPMI-1640

AL751G-20X50NO)

3. Washing RPMI 1640 media (1%FCS)

4. U-bottom 96 well culture plate

5. Ovalbumin (0µg/ml, 0.1µg/ml, 1.0µg/ml, 10µg/ml)

6. RBC lysis buffer

7. CO₂ Incubator

8. Frosted slides

2.6.3 Reverse Transcription and Polymerase Chain Reaction (RT-PCR):

For RNA isolation

1. Trizol or Tri reagent (cat# no. T9424-200ML, Sigma)

2. Diethyl pyrocarbon, DEPC (cat# no. D5758, Sigma)

3. Isopropanol

4. RNA precipitation solution

5. 75% Ethanol

6. Formamide
7. Chloroform
8. Cuvette for measuring O.D. at 260 nm and 280 nm
9. RNA isolation hood.

2.6.4 For Agarose gel electrophoresis:

1. RNA or DNA sample
2. 5% (for RNA) and 1% (for DNA, PCR products) Agarose gel
3. Autoclaved Mili Q water
4. 1X TAE buffer
4. Electrophoresis unit
5. Gene ruler 100bp DNA ladder, Fermentas (cat.# SM0241)
6. 6X loading dye
7. Gel Documentation system.
8. Gel loading dye (Bromo phenol blue)
9. Ethidium bromide (EtBr)

2.6.5 For cDNA synthesis:

1. RNA template
2. First Strand cDNA Synthesis Kit, Fermentas (cat. # K1612)
3. PCR machine (eppendorf thermal cycler)

2.6.6 For PCR amplification:

1. cDNA template
2. Forward IL-15 primer (sigma)
3. Reverse IL-15 primer (sigma)
4. Forward and reverse β -actin primer (eurofins mwg operon)

5. PCR Master Mix 2X, Fermentas (cat. # K0171)
6. PCR machine (Eppendorf)
7. Primer designing tool: IDT(Integrated DNA technology)

Web address: <http://eu.idtdna.com/analyzer/Applications/OligoAnalyzer>

2.7. Primers:

2.7.1. Mus musculus actin, beta (Actb), mRNA

NCBI Reference Sequence: NM_007393.3

>gi|145966868|ref|NM_007393.3| Mus musculus actin, beta (Actb), mRNA

Table 4: β -actin primer detail

Primer 5'→3'	Primer Sequence 5'→3'	%GC	T _m (°C)	Primer Length
Forward	TGGAATCCTGTGGCATCCATGAAAC	48	64	25
Reverse	TAAAACGCAGCTCAGTAACAGTCCG	48	64	25
Product Length: 349				

2.8.2. Mus musculus interleukin 15 (Il15), transcript variant 1, mRNA

NCBI Reference Sequence: NM_008357.2

Table 5: IL-15 primer detail

Primer 5'→3'	Sequence (5'→3')	Primer Start Position	%GC	T _M (°C)	Primer Length (Nucleotide)
Forward	GACAGTGACTTTCATCCCAG	732	50	61.1	20
Reverse	CCAAGTGGCTCATTATCTCC	1128	50	61	20
Product Size: 397 bp					

2.8 PROTOCOLS

2.8.1. RNA Isolation Protocol:

1. Preparation of tissue sample for RNA isolation:

- **100 mg** tissue was homogenized for 20-30 sec with **1 ml Tri reagent** in homogenizer.
- Centrifuged at 10,000 rpm for 10 min at 4°C
- Tri reagent is a denaturant, it denature the cell wall protein and lyse the cell.

2. Supernatant was taken and **250µl chloroform** was added. Mix well and incubate for 10 minutes at room temperature (RT).

- This help dissociating proteins present in cell extract.
- Centrifuge at 12,000 rpm for 15 min at 4°C

3. After centrifugation three phases will obtain. Take out the uppermost layer.

- The **upper aqueous layer** contains RNA, middle layer will contain DNA and lowest layer will contain protein.

Note: DNA and protein are extracted out in the organic phase while RNA remains in the aqueous phase.

4. Precipitation of RNA

- To the upper aqueous layer **200 µl RNA precipitating solution + 250 µl Isopropanol** was added.
- Mix gently for 2-3 minutes for 10 min at RT
- Centrifuged at 10,000 rpm for 15 min at 4°C
- This would precipitate out **RNA in the pellet form.**

4. The pellet thus obtained was washed twice with 1 ml of **75% ethanol** (prepared in DEPC treated water) by centrifuging it at 8000 rpm for 5 min at 4°C

- This wash is given to remove excess Isopropanol used for RNA precipitation.

5. Supernatant was discarded pellet obtained was allowed to air dry.

6. Finally pellet was dissolved in **20 µl formamide** and stored at -20°C or in liquid nitrogen.

Facts:

1. The yield of total RNA depends on the tissue or cell source, but it is generally 4-7 µg/mg starting tissue or 5-10 µg/10⁶ cells.

2. Dissolve the precipitated RNA in deionized formamide and store at -20°C. Formamide provides a chemically stable environment that also protects RNA against degradation by RNases. Purified, salt free RNA dissolves quickly in formamide

up to concentrations, samples of the RNA can be analyzed directly by gel electrophoresis, RT-PCR, or RNase protection saving time and avoiding potential degradation. (Chomczynski,1992)

3. If necessary, RNA can be recovered from formamide by precipitation with 4 volumes of ethanol or by diluting the formamide fourfold with 0.2N NaCl & then adding the 2 volumes of ethanol.

2.8.1.1 Agarose Gel Electrophoresis:

1.5% agarose gel is suitable for resolving RNA.

- 1.5% agarose gel was prepared in 1X TAE.
- 1X TAE was used as electrophoresis buffer
- Electrophoresis unit was installed and loaded with 1X TAE
- Sample was loaded with Bromophenol Blue dye in the proportion of 2:1
- Gel was run at 100V 20 mins.
- Gel picture were taken in Gel documentation system.

2.8.1.2 Spectrophotometric analysis:

The absorbance of RNA sample was measured at 260 nm and 280 nm. The nucleic acid concentration was calculated using the Beer-Lambert law, which predicts a linear change in absorbance with concentration.

$$A = \epsilon Cl,$$

Where, A= absorbance at a particular wavelength

C= concentration of nucleic acid

L= length of the spectrophotometric cuvette (1 cm)

E= extinction coefficient (ϵ for RNA is $0.025(\text{mg/ml})^{-1}\text{cm}^{-1}$)

For RNA Absorbance of 1 at 260 nm (A_{260} value 1) corresponds to $40\mu\text{g/ml}$ of single stranded RNA. A_{280} an idea about the amount of protein present in isolated RNA. The A_{260}/A_{280} ratio is used to assess RNA purity. An A_{260}/A_{280} ratio of 1.82 is indicative of highly purified RNA. RNA samples were diluted with DEPC- water and A_{260} was measured. And concentration was calculated using following formula:

$$\text{RNA concentration } (\mu\text{g}/\mu\text{l}) = \frac{A_{260} \times \text{Dilution Factor} \times 40}{1000}$$

1000

random hexamer are the primers that can be used for cDNA synthesis. Oligo dT primers specifically binds to poly-A tailed mRNA, while random hexamer non-specifically binds to all those RNA which are having sequence complementary to random hexamer.

Composition of Reaction 2.8.2 cDNA Synthesis:

RT-PCR (reverse transcription-polymerase chain reaction) is a technique used to study the expression of gene at mRNA levels. It is the most sensitive technique for mRNA detection and quantitation.

Reverse transcriptase is a RNA-dependent DNA polymerase enzyme isolated from retrovirus. Oligo dT and **Mixture for RT**

For cDNA synthesis various components were added following order

Table 6: Composition of Reaction Mixture for RT

S.no.	Reagents	Volume
1	Template DNA	*
2	Oligo dT primer	1 μ l
3	nuclease-free Water	Upto 11 μ l
	Total	11 μ l

*Volume of template RNA is adjusted according to RNA concentration to be used for cDNA synthesis

	Reagents	Volume
4	5X reaction Buffer	4 μ l
5	Ribolock RNase Inhibitor	1 μ l
6	10mM dNTP Mix	2 μ l
7	M-MuLV Reverse Transcriptase	2 μ l
	total volume	20 μ l

For cDNA synthesis tubes containing above mixture were incubated for 60min at 37⁰C (If oligo dT is used) and reaction was terminated by incubating at 70⁰C for 5min.

If random hexamer is used above mixture should be first incubated at 25⁰C for 5min followed by 60 min incubation at 37⁰C

2.8.3 Polymerase chain reaction (PCR):

Composition of Reaction Mixture for PCR

Table 7: Composition of Reaction Mixture for PCR

S.no.	Reagent	volume
1	PCR master mix	12.5
2	Template cDNA	2µl
3	Forward primer	1 µl
4	Reverse Primer	1 µl
5	Nuclease free water	Upto 25µl

After mixing all reagents give a short spin for proper mixing of all components. Above tube was placed in thermal cycler having with below mentioned PCR profile.

PCR profile for β-actin

Table 8: PCR profile for β-actin

S no.	Step	Temperature (°C)	Time	Number of cycle
1	Initial denaturation	95	3 min	1
2	Denaturation	95	30 sec	35
3	Annealing	56	45 sec	
4	Extension	72	60 sec	
5	Final extension	72	10 min	1

PCR profile for IL-15

Table 9: PCR profile for IL-15

S no.	Step	Temperature °C)	Time	1
1	Initial denaturation	95	3 min	35
2	Denaturation	55.8	30 sec	

3	Annealing	72	45 sec	1
4	Extension	72	60 sec	
5	Final extension	95	10 min	1

2.8.3 IL-15 induction:

In vitro induction of IL-15:

1. Isolation of spleenocytes was done by mincing through two frosted slides with 1X PBS and 1% glucose in petriplates
- Cells were spin down at 1000rpm for 10 mins, pellet was obtained.
2. Supernatant was discard, **pellet** was dissolved in 3 ml complete **RPMI media**.
 3. Cells were divided in required number of wells, and RPMI and **10% FCS** were added in each well.
- Note: FCS promotes and increases the degree of adherence of cells onto culture plate.
4. Cells were kept for incubation for **4 hours** at 37°C and 5% CO₂
- Note: Lymphocytes would be in the suspension while macrophage and DC adhered to culture plate.
5. **Lymphocyte suspension** was carefully removed and extra RPMI was added for culturing of Macrophage and DC.
 6. Cells were incubated **overnight** at 37°C, 5% CO₂.
- Note: Macrophage will be adhered and DC in the suspension.
7. **DC were removed** from the suspension by repeating a wash with PBS and transferred to sterile tube.
 8. To the Macrophage was added 0.5ml **Trypsin-EDTA** per well and kept for **3** minutes.
 9. Reaction was stopped by adding 3 ml RPMI.
 10. Macrophage were transferred into sterile tube and centrifuged at 1000 rpm for 10 mins to pellet down the cells.
 12. Pellet was dissolved in minimum volume of RPMI and taken the cell for counting.

Spleen was dissected cells were harvested and total RNA was isolated¹³. According to the cell count, pellet was dissolved in desired volume of media. Then the cells were divided in the required number of tubes.

14. Cells were stimulated with **LPS-10ng/ml** and **Poly I:C- 50ug/ml**.

15. Cells were incubated for **4 hours** at 37°C and 5% CO₂.

16. Cells were obtained in pellet by centrifuging at 1000rpm for 10mins.

17. Supernatant was discarded and **pellet** was dissolved in left volume of supernatant.

18. Adding 1 ml **Trizol** to the tube and cells were transferred in the eppendorf and proceeded for **RNA isolation** by Trizol method.

2.8.4 Enzyme-Linked Immunosorbent Assay (ELISA):

ELISA is a powerful and sensitive technique for the detection of antigen or antibody. The principle of this is same as that of radio immuno assay (RIA) i.e. an enzyme conjugated with an antibody reacts with a colourless substrate to generate a coloured reaction product and intensity of colour thus produced will give the concentration of antigen in the sample.

A number of variations of ELISA have been developed, allowing qualitative detection or quantitative measurement of either antigen or antibody, like direct ELISA, indirect ELISA, sandwich ELISA.

At a time interval of 24 hours and 72 hours after the immunization ~0.5 ml of blood was collected from three mice of each group by retro-orbital bleeding. Blood was allowed to clot for 4 hours at room temperature and then incubated at 4⁰C for better recovery of serum.

Serum was collected after centrifugation at 1500 rpm for 7 min. and IL-15 ELISA was performed using Ray biotech IL-15 ELISA kit.

RESULTS AND DISCUSSION

In the natural condition host might be encountered by more than one organism, so in this study we want to check the influence of one antigen on the immune response generated against the other antigen. PPD (purified protein derivatives) is a major mycobacterial protein antigen of CFA. Not much study is done about the role of PPD in the adjuvant activity of CFA. So here our objective is to study how PPD influencing the immune response generated against Ova.

Multivalent (polyvalent vaccine) vaccine is designed to elicit an immune response either to more than one infectious agent or to several different antigenic determinants of a single antigen. In natural condition when an organism is infecting any host, host immune system will be encountered with different antigens of that organism. So our next aim was to study whether two antigenic protein from same origin influence immune response generated by one antigenic protein. For that we used Ova (chicken Ovalbumin) and HEL (hen egg lysozyme) from the same source i.e. egg.

- **PPD shown enhanced antigen (Ova) specific CD4⁺ and CD8⁺T-cell responses.**

Eleven days after the post immunization cells were harvested from draining lymph node and were labelled with CFSE dye and T-cells were challenged *in vitro* with increasing concentration of Ova protein (0, 0.1,1.0 ,10 µg/ml) . Following 72 hr stimulation, CD4⁺ and CD8⁺ T-cells were labelled with anti-CD4 APC and anti-CD8 PE and analysed by Flow Cytometry.

Scheme of analysis:








Figure 11: Scheme of analysis



Ova concentration: 0µg/ml



Ova concentration: 0.1µg/ml



Ova concentration: 1.0µg/ml



Ova concentration: 10µg/ml









Figure 12: Mice were immunized with Ova, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0, 10 µg/ml)



Ova concentration: 0 μ g/ml



Ova concentration: 0.1 μ g/ml



Ova concentration: 1.0 μ g/ml



Ova concentration: 10 μ g/ml





Figure 13: Mice were immunized with Ova+PPD, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0, 10 μ g/ml)








Figure14: Mice were immunized with Ova+LPS, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0,10 µg/ml)



Ova concentration: 0µg/ml



Ova concentration: 0.1µg/ml



Ova concentration: 1.0µg/ml



Ova concentration: 10µg/ml





Figure 15: Mice were immunized with Ova+LPS+PPD, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0, 10 µg/ml)








Figure 16: Mice were immunized with Ova+ Zymosan, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0,10 µg/ml)


Ova concentration: 0 μ g/ml



Ova concentration: 0.1 μ g/ml



Ova concentration: 1.0 μ g/ml



Ova concentration: 10 μ g/ml










Figure 17: Mice were immunized with Ova+Zymosan +PPD, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0, 10 μ g/ml)



Ova concentration: 0 μ g/ml



Ova concentration: 0.1 μ g/ml



Ova concentration: 1.0 μ g/ml



Ova concentration: 10 μ g/ml





Figure 18: Mice were immunized with Ova+LPS+Zymosan, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0, 10 μ g/ml)




Figure 19: Mice were immunized with Ova+ LPS+ Zymosan+ PPD, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0,10 µg/ml)

PPD (purified protein derivatives) is a major constituent of CFA (Freund's complete antigen), CFA also contain TLR ligands like lipoarabinomannan (TLR2), heat shock protein (TLR4) etc. So we used zymosan (TLR2 ligand) and LPS (TLR4 ligand) for our study in different combination. Aim of our study was to show that how proteins of different origin can influence the immune response generated against another antigenic protein and whether PPD modulates TLR induce immune response or not. So we immunized mice with Ova as model antigen and measured CD4 T-cell proliferation by CFSE dilution assay.


Table: % of dividing CD4⁺ T-cells present in each division with respect to total dividing cell (TDC) and Total cell (dividing + Non dividing cells) and number of division cells has undergone upon *In vitro* challenge with Ova protein

A=% of dividing cells relative to total dividing cells; B=% of dividing cells relative Total cells (dividing + Non-dividing cells)

Combinations	Number of division cells has undergone	<i>In vitro</i> challenge with Ova (µg/ml)							
		0 µg/ml		0.1 µg/ml		1.0 µg/ml		10 µg/ml	
		A	B	A	B	A	B	A	B
Ova	1	75.0	0.9	87.3	2.6	91.7	4.5	3.9	100.0
	2	25.0	0.3	12.7	0.4	8.3	0.4	0.0	0.0
	3	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Ova+PPD	1	100.0	0.7	100.0	2.5	83.1	3.7	46.0	3.5
	2	0.0	0.0	0.0	0.0	16.9	0.8	29.9	2.3
	3	0.0	0.0	0.0	0.0	0.0	0.0	17.2	1.3
Ova +LPS	1	75.0	0.9	100.0	1.3	100.0	4.6	100.0	3.9
	2	25.0	0.3	0.0	0.0	0.0	0.0	0.0	0.0
	3	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Ova+LPS+PPD	1	100.0	0.9	100.0	1.8	100.0	2.0	100.0	8.3
	2	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Ova+Zymosan	1	100.0	2.5	86.4	3.0	100.0	4.4	100.0	3.6

	2	0.0	0.0	13.6	0.5	0.0	0.0	0.0	0.0
Ova+Zym+PPD	1	81.0	1.6	89.7	3.7	86.4	3.4	49.6	1.5
	2	19.0	0.4	10.3	0.4	13.6	0.5	28.1	0.9
	3	0.0	0.0	0.0	0.0	0.0	0.0	22.3	0.7
Ova+LPS+Zym	1	100.0	0.5	87.9	3.0	100.0	2.6	100.0	4.0
	2	0.0	0.0	12.1	0.4	0.0	0.0	0.0	0.0
Ova+LPS+Zym+PPD	1	100.0	1.4	100.0	2.5	83.9	3.7	100.0	5.5
	2	0.0	0.0	0.0	0.0	16.1	0.7	0.0	0.0

Figure 20: % of dividing Ova specific CD4⁺ T-cell after *in vitro* challenge with Ova ($\mu\text{g/ml}$): Mice were immunized with below mentioned combination and T-cell response was measured by CFSE dilution assay using flow cytometer




Values are mean \pm SD of 3 separate experiments (n=3mice per group)

ANOVA followed by Student's- Newman Keul's Test

*p<0.05 – As compared to Ova

§ p<0.05 – As compared to HEL



Values are mean± SD of 3 separate experiments (n=3mice per group)

ANOVA followed by Student's- Newman Keul's Test

@ p<0.05 – As compared to OZ

Figure 21: % of dividing Ova specific CD4⁺ T-cell after *in vitro* challenge with Ova (µg/ml) mice were immunized with below mentioned combination and T-cell response was measured by CFSE dilution assay using flow cytometer


Table: % of dividing CD8⁺ T- cells present in each division with respect to total dividing cell (TDC) and Total cell (dividing+ non dividing cells) and number of division cells has undergone upon *In vitro* challenge with Ova protein

		<i>In vitro</i> challenge with Ova (µg/ml)							
Combinations	Number of division cells has undergone	0 µg/ml		0.1 µg/ml		1.0 µg/ml		10 µg/ml	
		A	B	A	B	A	B	A	B
Ova	1	100.00	6.65	94.49	10.32	86.29	10.05	85.86	17.11
	2	0.00	0.00	5.51	0.60	13.71	1.60	10.77	2.15
Ova+PPD	1	58.61	10.40	94.34	17.09	86.97	23.62	83.82	23.62
	2	35.71	6.34	4.88	0.88	11.15	3.03	11.95	3.03
	3	4.78	0.85	0.78	0.14	1.88	0.51	4.23	0.51

	4	0.91	0.16	0.00	0.00	0.00	0.00	0.00	0.00
Ova +LPS	1	100.00	12.87	100.00	18.39	100.00	15.24	75.70	10.53
	2	0.00	0.00	0.00	0.00	0.00	0.00	24.30	3.38
Ova+LPS+PPD	1	100.00	3.14	70.41	2.73	90.45	4.25	90.05	3.21
	2	0.00	0.00	29.59	1.15	9.55	0.45	9.95	0.35
Ova +Zymosan	1	94.24	13.78	60.69	8.86	86.30	17.65	76.11	9.62
	2	5.76	0.84	30.34	4.43	12.09	2.47	22.67	2.86
	3	0.00	0.00	8.97	1.31	1.62	0.33	1.21	0.15
Ova+Zym+PPD	1	84.19	11.30	52.30	10.67	71.00	15.73	75.20	18.50
	2	15.81	2.12	42.31	8.63	27.67	6.13	24.80	6.10
	3	0.00	0.00	5.39	1.10	1.32	0.29	0.00	0.00
Ova+LPS+Zym osan	1	82.28	2.54	45.67	7.38	74.80	10.95	91.02	19.89
	2	17.72	0.55	46.72	7.55	17.57	2.57	7.68	1.68
	3	0.00	0.00	7.61	1.23	7.63	1.12	1.30	0.29
Ova+LPS+Zym +PPD	1	58.80	6.30	68.86	12.82	85.08	18.81	70.56	12.05
	2	33.63	3.60	26.80	4.99	12.60	2.79	26.94	4.60
	3	7.57	0.81	4.35	0.81	1.52	0.34	2.50	0.43
	4	0.00	0.00	0.00	0.00	0.80	0.18	0.00	0.00

A=% of dividing cells relative to total dividing cells; B==% of dividing cells relative to total cells (dividing+ Non-dividing cells)

Ova=Ovalbumin; LPS= Lipopolysaccharide; Zym=Zymosan; PPD= purified protein derivatives




Values are mean± SD of 3 separate experiments (n=3mice per group)

ANOVA followed by Student's- Newman Keul's Test

*p<0.05 – As compared to Ova

& p<0.05 – As compared to OL

Figure 22: % of dividing Ova specific CD8⁺ T-cell after *in vitro* challenge with Ova (µg/ml) mice were immunized with below mentioned combination and T-cell response was measured by CFSE dilution assay using flow cytometer



Values are mean± SD of 3 separate experiments (n=3mice per group)

ANOVA followed by Student's- Newman Keul's Test

p<0.05 – As compared to Ova

@ p<0.05 – As compared to OZ

Figure 23: % of dividing Ova specific CD8⁺ T-cell after *in vitro* challenge with Ova (µg/ml), Mice were immunized with below mentioned combination and T-cell response was measured by CFSE dilution assay using flow cytometer

- **PPD is enhancing T-cell proliferation in the presence and absence of TLR ligands**




Figure 24: % of dividing ova specific CD4⁺ T-cells at different *in vitro* Ova challenge.




Figure 25: % of dividing ova specific CD4⁺ T-cells at different *in vitro* Ova challenge.

From our study we found that PPD is positively modulating response generated against Ova. When Ova alone was injected cells has undergone only two division cycle and the % of responding cells in second cycle was only 25%, but when injected with PPD, presence of PPD has pushed cells 17% of responding cells to third division.

Similar is the case when Ova and LPS were injected with or without PPD. In the absence of PPD cells has undergone only two division cycle and only 13% of the responding cells were in second division, while in the presence of PPD has pushed the cells to 3 division cycle and the % of responding cells in third division was 23.3%

So our results suggest that PPD is enhancing TLR induced immune response and our study also shows that not only TLR ligands of CFA is responsible for adjuvant activity of CFA but PPD (purified protein derivatives) of *Mycobacterium tuberculosis* contributes majorly to make it as a strong adjuvant.

We further suggest to study molecular details that play a crucial role in modulating TLR induced response via PPD. It is known that all TLR ligands except TLR3 mediates their action via TRAF6 and NF-kB protein in MyD88 dependent and MyD88 independent pathway. So we assume that PPD modulates TLR induce response via this pathway and for that we suggest to measure TRAF6 protein level by immunoprecipitation.

- **PPD upregulates interleukin-15 (IL-15) expression, in the presence or absence of TLR ligand**

S. no.	Combinations	Concentration of IL-15 (pg/ml)	
		After 24h	After 72h
1	Ova	60.98	-
2	Ova+PPD	90.79	-
3	Ova+LPS+Zymosan	115.18	102.53
4	Ova+LPS+Zym+PPD	256.94	102.08





Figure 26: Amount of IL-15 produced after 24 of immunization, blood was collected 24 hours post immunization and serum was separated and IL-15 level was measured by using ELISA kit (Ray Biotech)

When only Ova was injected the amount of IL-15 produced was quite less as compared when Ova was injected with PPD. Similar is the results when Ova, LPS and Zymosan were injected with or without PPD, with PPD the amount of IL-15 produced was significantly higher than without PPD. This data suggest that PPD is upregulating IL-15 expression.


- **Hen egg lysozyme (HEL) injection modulates the CD4⁺ and CD8⁺ T-cell response generated against another antigen (Ova).**

Eleven days after the post immunization cells were harvested from draining lymph node and were labelled with CFSE dye and T-cells were challenged *in vitro* with increasing concentration of Ova protein (0, 0.1, 1.0, 10 µg/ml). Following 72 hr stimulation, CD4⁺ T-cells were labelled with anti-CD4 APC and analysed by Flow Cytometry


In vitro Ova concentration: 0 μ g/ml



In vitro Ova concentration: 0.1 μ g/ml



In vitro Ova concentration: 1.0 μ g/ml



In vitro Ova concentration: 10 μ g/ml




Figure 27: Mice were immunized with Ova, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0, 10 μ g/ml)




Figure 28: Mice were immunized with HEL, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0, 10 µg/ml)




Figure 29: Mice were immunized with Ova+ HEL, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0,10 µg/ml)




Figure 30: Mice were immunized with Ova+LPS+Zymosan, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0, 10 µg/ml)




Figure 31: Mice were immunized with Ova+HEL+LPS+ Zymosan, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0,10 µg/ml)




Figure 32: Mice were immunized with HEL+LPS+Zymosan, T-cell isolated from draining lymph node and challenged *in vitro* with Ova (0, 0.1, 1.0, 10µg/ml)


Multivalent (polyvalent vaccine) vaccine is designed to elicit an immune response either to more than one infectious agent or to several different antigenic determinants of a single antigen. So our next aim was to study whether two antigenic protein from same origin influence immune response generated by one antigenic protein. For that we used Ova (chicken Ovalbumin) and HEL (hen egg lysozyme) and immunize mice intraperitoneally with two TLR ligands (LPS and zymosan) and two antigenic proteins i.e. Ova and HEL and measured CD4 T-cell proliferation by CFSE dilution assay using flow cytometer.

Table: % of dividing CD8⁺ T-Cells present in each division with respect to total dividing cell (TDC) and Total cell (Dividing+ Non dividing cells) and number of division cells has undergone upon *In vitro* challenge with Ova protein.

A=% of dividing cells relative to total dividing cells; B==% of dividing cells relative to total cells (dividing+ Non-dividing cells)

Combination injected	Number of division	<i>In vitro</i> challenge with Ova (µg/ml)							
		0 µg/ml		0.1 µg/ml		1.0 µg/ml		10 µg/ml	
		A	B	A	B	A	B	A	B
HEL	1	77.42	1.82	95.75	15.47	75.68	2.67	78.37	5.82
	2	12.90	0.30	3.09	0.50	16.22	0.57	19.59	1.45
	3	9.68	0.23	0.77	0.12	5.41	0.19	1.63	0.12
	4	0.00	0.00	0.39	0.06	2.70	0.10	0.41	0.03
	5	0.00	0.00	0.00	0.00	0.23	0.01	0.00	0.00
	6	0.00	0.00	0.00	0.00	0.11	0.00	0.00	0.00
Ova+HEL	1	85.71	4.24	58.92	3.71	81.93	2.50	82.61	3.56
	2	10.71	0.53	27.62	1.74	16.87	0.52	17.39	0.75
	3	3.57	0.18	11.51	0.72	1.20	0.04	0.00	0.00
	4	0.00	0.00	1.96	0.12	0.00	0.00	0.00	0.00

HLZ	1	96.82	8.23	88.89	4.30	95.49	14.52	43.73	5.71
	2	3.18	0.27	11.11	0.54	4.14	0.63	56.27	7.34
	3	0.00	0.00	0.00	0.00	0.37	0.06	0.00	0.00
Ova+HEL+LPS+ Zymosan	1	75.60	8.47	74.63	7.36	75.80	5.20	78.20	9.28
	2	22.42	2.51	23.13	2.28	23.35	3.10	18.05	2.14
	3	1.54	0.17	2.24	0.22	0.85	0.20	2.79	0.33
	4	0.44	0.05	0.00	0.00	0.00	0.00	0.97	0.11
Ova+LPS+ Zymosan	1	98.98	9.22	100.00	13.73	98.67	11.62	100.00	3.47
	2	1.02	0.10	0.00	0.00	1.33	0.16	0.00	0.00



Values are mean \pm SD of 3 separate experiments (n=3 mice per group)

ANOVA followed by Student's- Newman Keul's Test

*p<0.05 – As compared to Ova

\$ p<0.05 – As compared to HEL

Figure 33: % of dividing Ova specific CD8⁺ T-cell after *in vitro* challenge with Ova (µg/ml), mice were immunized with below mentioned combination and T-cell response was measured by CFSE dilution assay using flow cytometer


Table: % of dividing CD4⁺ T-Cells present in each division with respect to total dividing cell (TDC) and Total cell (dividing+ non dividing cells) and number of division cells has undergone upon *In vitro* challenge with Ova protein.

A=% of dividing cells relative to total dividing cells; B==% of dividing cells relative to total cells (dividing+ Non-dividing cells)

Combination	Number of division	<i>In vitro</i> challenge with Ova ($\mu\text{g/ml}$)							
		0 $\mu\text{g/ml}$		0.1 $\mu\text{g/ml}$		1.0 $\mu\text{g/ml}$		10 $\mu\text{g/ml}$	
		A	B	A	B	A	B	A	B
Ova	1	100	2.07	100	12.79	87.17	7.9	100	3.18
	2	0.00	0.00	0.00	0.00	12.82	4.23	0.00	0.00
	3	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
HEL	1	79.75	13.79	56.44	14.31	52.66	11.64	29.22	7.6
	2	14.56	2.52	22.09	5.60	25.08	5.54	58.90	15.67
	3	4.75	0.82	17.18	4.35	16.93	3.74	9.36	2.4
	4	0.95	0.16	4.29	1.09	5.33	1.18	2.51	0.66
	5	0.00	0.00	0.46	0.12	2.12	0.47	0.86	0.22
	6	0.00	0.00	0.31	0.08	0.18	0.04	0.23	0.06
Ova+HEL	1	56.54	8.07	41.92	7.42	54.61	6.10	100.00	13.25
	2	43.46	6.20	42.24	7.48	40.59	4.53	0.00	0.00
	3	0.00	0.00	15.84	2.80	4.80	0.54	0.00	0.00
	4	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
HLZ	1	61.73	18.38	100.00	21.59	100.00	13.89	100.00	20.03
	2	38.27	11.39	0.00	0.00	0.00	0.00	0.00	0.00
	3	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
Ova+HEL+LPS+	1	64.22	10.79	39.77	4.42	40.3	3.1	65.16	8.86

Zymosan	2	19.85	3.34	39.02	4.34	42.5	3.6	18.36	2.50
	3	13.73	2.31	18.01	2.00	13.2	1.2	12.93	1.76
	4	2.21	0.37	3.19	0.35	3.6	0.3	3.56	0.48
	5	0.00	0.00	1.69	0.19	0.4	0.2	2.68	0.36
	6	0.00	0.00	0.00	0.00	0.00	0.00	0.86	0.12
Ova+LPS+ Zymosan	1	65.54	7.72	62.81	11.07	55.85	8.92	64.03	9.53
	2	20.95	2.47	24.76	4.36	23.04	3.68	29.48	4.39
	3	9.46	1.11	8.71	1.53	15.56	2.48	5.07	0.75
	4	4.05	0.48	3.72	0.66	5.55	0.89	1.43	0.21
	5	2.45	0.29	0.52	0.09	0.88	0.14	0.00	0.00
	6	0.44	0.05	0.36	0.06	0.00	0.00	0.00	0.00

Figure 34: % of dividing Ova specific CD4⁺ T-cell after *in vitro* challenge with Ova ($\mu\text{g/ml}$), mice were immunized with below mentioned combination and T-cell response was measured by CFSE dilution assay using flow cytometer



Values are mean \pm SD of 3 separate experiments (n=3mice per group)

ANOVA followed by Student's- Newman Keul's Test

*p<0.05 – As compared to Ova

p<0.05 – As compared to OLZ

- ***In vitro* Ova challenge is helping more number of cells to undergo subsequent cell division cycle**




Figure 35: % of division Ova specific CD4⁺ T-cells has undergone at 0 µg/ml *in vitro* Ova concentration (each colour indicate one division cycle)




Figure 36: % of division Ova specific CD4⁺ T-cells has undergone at 1.0 µg/ml *in vitro* Ova concentration (each colour indicate one division cycle)

- **Increasing *in vitro* Ova challenge in presence of TLR ligands, decreases T-cell proliferation**




Figure 37: % of dividing Ova specific CD4⁺ T-cells with increase in *in vitro* Ova challenge

Most of the Ova specific T-cells taken from mice immunized with Ova alone get stuck in first division even upon the *in vitro* challenge with Ova protein, but in the HEL alone as *in vitro* challenge by the Ova protein increases in the concentration from 0 to 10 µg/ml, it is helping cells to undergo the further division cycles upto 6 cycle, which suggest that *in vitro* challenge by Ova might be upregulating the MHC molecule which present more HEL specific epitope to specific T-cell, or it is increasing the contact time between MHC-epitope complex and TCR (T-cell receptor) and leads to higher proliferation or it may be due to cross-reactivity of Ova protein. One another possibility may be that after *in vitro* challenge by Ova, it is engulfed and processed by APCs (antigen presenting cells) and providing such specific cytokines rich niche i.e. signal 3 which in turn leads to T-cell proliferation. But when both Ova and HEL together are injected *in vivo* we observed that now cells are not pushed much for proliferation like in case of HEL alone but it is limited to 3 division cycles only, which might be due to dilution of cytokines signals i.e. signal 3 which is generated after *in vitro* challenge by Ova protein, but at the same time it is sufficient enough for more proliferation compared to the Ova alone.

IL-15 ELISA:

S. no.	Combination injected	Concentration of IL-15 (pg/ml) At 72 h post immunization
1	HEL	95.53
2	Ova+HEL	87.40
3	HEL+LPS+Zym	171.64
4	Ova+HEL+LPS+Zym	71.14




Figure 38: Amount of IL-15 produced after 72 post immunization

The difference in the amount of IL-15 produced in HEL, OH and OLZ was not quite high, but the difference is very significant between HLZ and OHLZ. When Ova is injected with HLZ the IL-15 level has reduced significantly, suggesting that there are some unknown interaction between HEL and Ova.

Further our aim was to co-relate the proliferation of T-cell to innate immunity. The APCs (Antigen presenting cells) process antigen and present epitopes to naïve T-cells. During this phase many immune as well as non-immune cells produce cytokines such as IL-2, IL-7, IL-15 etc, which play an important role in T-cell proliferation and survival at many stages of immune response. To investigate the expression of IL-15 we isolated mRNA from splenic dendritic cell stimulated with TLR ligands and performed RT-PCR.

3.2 RT-PCR

3.2.1. Standardization of RNA isolation by Trizol method:

RNA was isolated from spleen sample by Trizol method. To check the integrity of isolated RNA electrophoresis was performed using 1.5% Agarose gel electrophoresis.




Figure 39: Electrophoresis of isolated RNA sample

Ideally 28S, 18S and 5S rRNA should be observed in the Agarose gel. And according to the standard observations, the intensity of 28S rRNA should be highest and that of 18S and 5S rRNA should in decreasing manner.

The expectation was matched in the above gel-electrophoresis, suggesting that RNA isolated by our protocol were of good quality. The presence of intact 28S and 18S rRNA band indicate that the mRNA in the total sample is intact and not degraded by RNase.

3.2.2 Spectrophotometric analysis of RNA sample:

The concentration and purity of RNA sample was determined spectrophotometrically taking A₂₆₀/A₂₈₀nm ratio. Nucleotides give maximum absorbance a 260 nm(A₂₆₀) it reflects the concentration of RNA in our sample, While protein gives maximum absorbance at 280nm (A₂₈₀nm). Ratio of absorbance at 260 and 280nm (A₂₆₀/A₂₈₀)

indicates the purity of RNA sample. A_{260}/A_{280} ratio in the range of 1.8-2.0 indicates a highly purified RNA preparation.

Absorbance and concentration of above mentioned RNA sample is given below:

Sr. No.	Sample	Absorbance at 260nm A_{260nm}	Absorbance at 280nm A_{280nm}	A_{260}/A_{280nm}	Concentration ($\mu\text{g/ml}$)
1.	Sample 1	0.126	0.078	1.62	2.52
2.	Sample 2	0.259	0.125	2.07	5.18
3.	Sample 3	0.217	0.129	1.68	4.34
4.	Sample 4	0.298	0.142	2.09	5.96
5.	Sample 5	0.281	0.137	2.05	5.62
6.	Sample 6	0.286	0.172	1.66	5.72

Dilution factor- 500

Table: Absorbance at 260nm, 280nm and A_{260}/A_{280} .

3.3 Standardization of cDNA synthesis:




Figure 40: electrophoresis of cDNA transcribed from isolated RNA

cDNA synthesis was carried out for the RNA sample at 37°C for 1 hour where we could successfully Reverse Transcribe the RNA sample of 1.33 µg/ml into cDNA and was observed by gel electrophoresis using 1.2% Agarose gel.

3.4 Standardization of PCR amplification:

PCR was carried out using **control gene, β actin**. Primers for β-actin were designed whose product size was **349bp**.




Figure 41: Electrophoresis of PCR for β-actin at 56°C

PCR for control gene (**β-actin**) was performed at **56°C** where we could successfully obtain amplified bands for RNA sample with concentration **0.20 ug** (lane 2 and 3) with template volume 1 ul and 5 ul. Electrophoresis was carried out using 1% agarose gel.

3.5 *In vitro* induction of IL-15:

Dendritic cells were isolated from spleen using protocol as given above. Dendritic cells were stimulated with **LPS 10 ng/ml** and **Poly I:C 50 ug/ml**, **LPS(5ng/ml) + Poly I:C (25ug/ml)**. After 4 hours RNA was isolated by the protocol as given above. The electrophoresis of the isolated RNA samples was carried out using 1.5% Agarose gel electrophoresis.

a. RNA isolation from splenic dendritic cells stimulated with TLR ligands




Figure 42: RNA from splenic dendritic cells

Ideally 28S rRNA and 18S rRNA band should be observed in the agarose gel. According to the standard observation, the intensity of 28S rRNA is higher than the 18S rRNA. The expectation is matched in the above agarose gel, where we could successfully isolate RNA from **control** (lane 1) and **LPS** (lane 2), **Poly I:C** (lane 3), **LPS+Poly I:C** (lane 4), **LPS II** (lane 5) induced sample.

.5.1 Polymerase Chain Reaction of DC cells stimulated with TLR ligands (*in vitro*)




Figure 43: PCR run for β-actin

PCR was carried out for control gene (**β -actin**) at 56°C. Electrophoresis was carried out using 1% agarose gel where we could successfully obtain band in both, **control** (lane 2) and **LPS** (lane 3) induced sample run against the **ladder** of **100bp**. Our control gene has the product size 349 bp, which we could match against the ladder.

SUMMARY

In the natural condition host might be encountered by more than one organism, so in this study we want to check the influence of one antigen on the immune response generated against the other antigen. PPD (purified protein derivatives) is a major mycobacterial protein antigen of CFA. Not much study is done about the role of PPD in the adjuvant activity of CFA. So here our objective is to study how PPD influencing the immune response generated against Ova. And to test this idea whether one antigen from pathogen can influence the immunogenicity of another antigen, we developed a model that mimic the infection and measured the T-cell response against the targeted antigen. We found that PPD enhances antigen (Ova) specific CD4⁺ and CD8⁺T-cell responses. And our study also shows that not only TLR ligands of CFA is responsible for adjuvant activity of CFA but PPD (purified protein derivatives) of *Mycobacterium tuberculosis* contributes majorly to make it as a strong adjuvant. We further suggest to work on the molecular details of the immunomodulatory activity of PPD to induce strong cellular responses

Multivalent (polyvalent vaccine) vaccine is designed to elicit an immune response either to more than one infectious agent or to several different antigenic determinants of a single antigen. In natural condition when an organism is infecting any host, host immune system will be encountered with different antigens of that organism. So our next aim was to study whether two antigenic protein from same origin influence immune response generated by one antigenic protein. For that we used Ova (chicken Ovalbumin) and HEL (hen egg lysozyme) from the same source i.e. egg. We found that Hen egg lysozyme (HEL) injection modulates the CD4⁺ and CD8⁺ T-cell response generated against another antigen (Ova). We also found that HEL alone is helping more number of the cells for further division compared to Ova alone, which might be due the *in vitro* challenge by Ova protein provides cytokines rich niche for T-cell proliferation or due to the cross-reactivity of Ova protein. But when both Ova and HEL injected *in vivo* cell proliferation is less in comparison to HEL alone.

Further our aim was to co-relate the proliferation of T-cell to innate immunity. We have standardized the protocol for RNA isolation and got high quality RNA. We have standardized the conditions of reverse transcriptase (RT)-PCR for control gene, β -actin. To standardize the condition for IL-15, we performed cDNA synthesis with different RNA concentration and with various template concentrations. We also measured systemic IL-15 level using ELISA-kit, collecting blood 24h and 72h post immunization. And we found that PPD upregulates interleukin-15 (IL-15) expression, in the presence or absence of TLR ligand.

REFERENCES

1. Akira S. et al. "Pathogen Recognition and Innate Immunity." *Cell* 124 (4): 783-801 (2006).
2. Alpdogan, Ö. and Brink M. "IL-7 and IL-15: therapeutic cytokines for immunodeficiency." *Trends in Immunology* 26 (1): 56-64, (2005).
3. Alves N. et al. "IL-15 induces antigen-independent expansion and differentiation of human naive CD8+ T cells *in vitro*." *Blood* 102 (7): 2541 -2546 (2003).
4. Abul K. "Cellular and Molecular Immunology, Saunders Elsevier, India, 189-191, (2007).
5. Berard M. et al. "IL-15 Promotes the Survival of Naive and Memory Phenotype CD8+ T Cells." *The Journal of Immunology* 170 (10): 5018 -5026 (2003).
6. Billiau A. and Matthys P. "Modes of action of Freund's adjuvants in experimental models of autoimmune diseases." *Journal of Leukocyte Biology* 70 (6): 849 -860 (2001).
7. Bulfone-paus S. et al. "Death deflected: IL-15 inhibits TNF- α -mediated apoptosis in fibroblasts by TRAF2 recruitment to the IL-15R α chain." *The FASEB Journal* 13 (12): 1575 -1585 (1999).
8. Byrd-Leifer, C.A. et al.. "The role of MyD88 and TLR4 in the LPS-mimetic activity of Taxol." *European Journal of Immunology* 31 (8): 2448-2457 (2001).
9. CASTRO AP. et al.. "Adjuvant treatment increases the resistance to Mycobacterium avium infection of Mycobacteria-susceptible BALB/c mice." *Clinical & Experimental Immunology* 92 (3): 466-472 (1993).
10. David F. Tough et al. "T Cell stimulation *in vivo* by Lipopolysaccharide (LPS)" . *JEM* 185 (12): 2089-2094 (1997).
11. De Silva BS et al. . "Purified Protein Derivative (PPD) as an Immunogen Carrier Elicits High Antigen Specificity to Haptens." *Bioconjugate Chemistry* 10 (3): 496-501 (1999) .
12. Fehniger TA and Caligiuri MA. "Interleukin 15: biology and relevance to human disease." *Blood* 97 (1): 14 -32 (2001).
13. Frasnelli M. et al. "TLR2 modulates inflammation in zymosan-induced arthritis in mice." *Arthritis Res Ther* 7 (2): R370-R379 (2005).

14. Galelli A, Chedid L. "Modulation of myelopoiesis in vivo by synthetic adjuvant-active muramyl peptides: induction of colony-stimulating activity and stimulation of stem cell proliferation." *Infect. Immun.* 42, 1081–1085 (1983).
15. Giri, J.G. et al. (1995). "Identification and cloning of a novel IL-15 binding protein that is structurally related to the α -chain of the IL-2 receptor." *EMBO J.* 14, 3654–3663.
16. Golenbock, D. T. and M. J. Fenton (2001). "Extolling the diversity of bacterial endotoxins." 2 (4): 286-288.
17. Goodwin, R. G. and D. Friend, et al. (1990). "Cloning of the human and murine interleukin-7 receptors: Demonstration of a soluble form and homology to a new receptor superfamily." *Cell* 60 (6): 941-951.
18. Grabstein, K. H. and J. Eisenman, et al. (1994). "Cloning of a T cell growth factor that interacts with the beta chain of the interleukin-2 receptor." *Science* 264 (5161): 965-968.
19. Herbert, W. J. (1968). "The mode of action of mineral-oil emulsion adjuvants on antibody production in mice." *Immunology* 14, 301–318.
20. Kaech, S. M. and S. Hemby, et al. (2002). "Molecular and Functional Profiling of Memory CD8+ T Cell Differentiation." *Cell* 111 (6): 837-851.
21. Kaech, S. M. and J. T. Tan, et al. (2003). "Selective expression of the interleukin 7 receptor identifies effector CD8+T cells that give rise to long-lived memory cells." 4 (12): 1191-1198.
22. Kol, A., Lichtman, A. H., Finberg, R. W., Libby, P., Kurt-Jones, E. A. (2000). "Heat shock protein (HSP) 60 activates the innate immune response: CD14 is an essential receptor for HSP60 activation of mononuclear cells." *J. Immunol.* 164, 13–17.
23. Maki, K. and S. Sunaga, et al. (1996). "Interleukin 7 receptor-deficient mice lack $\gamma\delta$ - T cells." *Proceedings of the National Academy of Sciences* 93 (14): 7172-7177.
24. Mattei, F. and G. Schiavoni, et al. (2001). "IL-15 Is Expressed by Dendritic Cells in Response to Type I IFN, Double-Stranded RNA, or Lipopolysaccharide and Promotes Dendritic Cell Activation." *The Journal of Immunology* 167 (3): 1179-1187.
25. Muegge, K. and M. P. Vila, et al. (1993). "Interleukin-7: a cofactor for V(D)J rearrangement of the T cell receptor beta gene." *Science* 261 (5117): 93-95.
26. MÜSSENER, Å. and L. KLARESKOG, et al. (1995). "TNF- α Dominates

- Cytokine mRNA Expression in Lymphoid Tissues of Rats Developing Collagen- and Oil-Induced Arthritis." *Scandinavian Journal of Immunology* 42 (1): 128-134.
27. Napolitano, L. A. and R. M. Grant, et al. (2001). "Increased production of IL-7 accompanies HIV-1-mediated T-cell depletion: implications for T-cell homeostasis." *J. Immunol.* 167 (1): 73-79.
28. Ohashi, K., Burkart, V., Flohe, S., Kolb, H. (2000). "Heat shock protein 60 is a putative endogenous ligand of the toll-like receptor-4 complex." *J. Immunol.* 164, 558–561.
29. Sawa, Y. and Y. Arima, et al. (2009). "Hepatic Interleukin-7 Expression Regulates T Cell Responses." *Immunity* 30 (3): 447-457
30. Saitoh, S. I. and S. Akashi, et al. (2004). "Lipid A antagonist, lipid IVa, is distinct from lipid A in interaction with Toll like receptor 4 (TLR4) MD 2 and Ligand induced TLR4 oligomerization." *International Immunology* 16 (7): 961 -969
31. Saitoh, S. and K. Miyake (2009). "Regulatory molecules required for nucleotide-sensing Toll-like receptors." *Immunological Reviews* 227 (1): 32-43.
32. Sato, M. and H. Sano, et al. (2003). "Direct Binding of Toll-Like Receptor 2 to Zymosan, and Zymosan-Induced NF- κ B Activation and TNF- α Secretion Are Down-Regulated by lung Collectin Surfactant Protein A." *The Journal of Immunology* 171 (1): 417 -425.
33. Schluns, K. S. and W. C. Kieper, et al. (2000). "Interleukin-7 mediates the homeostasis of naive and memory CD8+T cells in vivo." *J. Immunol.* 165 (5): 426-432.
34. Schluns, K. S. and L. Lefrancois (2003). "Cytokine control of memory T-cell development and survival." *J. Immunol.* 171 (4): 269-279.
35. Takeda, K. and S. Akira (2005). "Toll-like receptors in innate immunity." *International Immunology* 17 (1): 1 -14.
36. Varas, A. and A. Vicente, et al. (1998). "Interleukin-7 Influences the Development of Thymic Dendritic Cells." *Blood* 92 (1): 93 -100.
37. Von Freeden-Jeffry, U. and N. Solvason, et al. (1997). "The Earliest T lineage Committed Cells Depend on IL-7 for Bcl-2 Expression and Normal Cell Cycle Progression." *Immunity* 7 (1): 147-154.
38. Weiner, G. J. (2000). "The immunobiology and clinical potential of immunostimulatory CpG oligodeoxynucleotides." *Journal of Leukocyte Biology* 68 (4): 455 -463.
39. Wilfinger WW, Mackey K, and Chomczynski P (1997) Effect of pH and ionic strength on the spectrophotometric assessment.

APPENDIX

- **DEPC- treated water:**

Prepare 0.1% v/v DEPC added in to sterile Nuclease Free Water

(Milipore). Keep for 1 hour at 37°C until the dissolution of DEPC in water.

- **RNA Precipitation Solution:**

1.2 M NaCl

0.8M Trisodium Citrate

Dissolve in DEPC-treated water.

- **50 X TAE:**

Tris acetate-EDTA: (PH: 8.4)

Tris base: 24.20 g

Glacial acetic acid: 5.7 ml

EDTA: 3.72 g

Prepared in sterile Nuclease Free Water (Milipore).

- **Gel loading dye:**

Preparation for 10 ml dye

50% glycerol -5 ml

10mM EDTA(ph-8.0)- 37mg

0.25%w/v BPB- 25 mg

- **Ethidium Bromide:**

EtBr: 2mg/ml

Prepare in DEPC-treated water, protect from light and store at 4°C

- **Phosphate Buffered Saline(PBS): (1X)**

PH: 7.4

Preparation for 300 ml

NaCl- 2.4 g

KCl- 60 mg

Na₂HPO₄: 432 mg

KH₂PO₄: 72 mg

- **RBC lysis buffer:**

Preparation for 200ml

KHCO₃: 0.2 g

NH₄Cl: 1.658 g

EDTA: 7.4 mg

- **CFSE:**

Stock: 5mM

Stock was prepared by dissolving 2.78 mg CFSE in 1 ml Dimethyl Sulphoxide (DMSO).

- **PHA:**

Stock: 1mg/ml

Stock was prepared by dissolving 1 mg in 1 ml sterile Nuclease Free water (milipore)